


SEE SHEET 2 OF 2 FOR COMPONENT DETAILS

CAPTIVE 4-40 CAP SCREWS FOR SECURING
 $\phi 1"$ OAP TO MOUNTING PLATE
 2.0mm HEX 3 PLACES




$\phi 6.0\text{mm}$ AGE ROD BORE
 6.4mm [0.25in] DEEP
 FOR USE WITH ER-SERIES RODS
 (NOT INCLUDED)
 4 PLACES, BOTH SIDES


M4 X 0.7 LOCKING SCREW
 2.0mm [5/64in] HEX
 4 PLACES, BOTH SIDES

ARROW ENGRAVING INDICATES
 COLLIMATED BEAM AXIS

M6 x 1.0 MOUNTING HOLE
 7.6mm [0.30in] DEEP
 BOTH SIDES




J11 (1.035-40) SERIES
 INTERNAL THREAD
 5.8mm [0.23in] DEEP
 $\phi 25.4\text{mm}$ [1.00in]
 CLEAR APERTURE



LOCKING SCREW FOR SECURING
 OPTIC PLATE TO CORNER CUBE
 2.0mm [5/64in] HEX








LASER ENGRAVED BOX FOR
 LABELING OF MOUNTED OPTIC



3/16-100
 FINE ADJUSTMENT SCREW
 WITH REMOVABLE KNOB
 2.0mm [5/64in] HEX
 3 PLACES

FOR INFORMATION ONLY
 NOT FOR MANUFACTURING PURPOSES

SEE SHEET 1 OF 2 FOR ASSEMBLY DETAILS



FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES